


HFRR Humidity Cabinet; Temperature Control Unit (TCU)

Changes to the ISO 12156 Standard

The TCU has been developed to ensure all our customers can meet the revised ISO 12156 requirements which are due for release later this year. The changes to the standard include **tighter** limits for temperature and humidity ranges for the ambient lab conditions, as shown in Figure 1. The hatched area displays the lab conditions in which the TCU can operate effectively and reach the required testing conditions in a suitable amount of time.

To follow the progress of the ISO revision please see the ISO website: <https://www.iso.org/standard/72481.html>





Purchasing options

Below is a comprehensive list of all the options we offer for the humidity cabinet and TCU. For further information and pricing please contact your local distributor.

Part Number	Title
PCS-003512	Complete Humidity Cabinet with TCU
PCS-003534	Retrofit Temperature Control Unit
PCS-002799 (HFRHCAB)	Humidity controlled cabinet for HFRR
PCS-003553	Temperature Control Unit