


Microscope Holder for MTM Specimens



The microscope holder brings the MTM specimens into the correct angle for an effortless focus on the contact track / wear scar. The specimen on the holder can be rotated easily by hand.

Available for both the MTM1 and MTM2

Supplied Parts

Mechanical - Holder, Barrel Adapter and Securing Screw

Software – None

Microscope Holder for MTM Specimens



Microscope kit + a MTM steel barrel and adapter



$\frac{3}{4}$ " ball on holder and $\frac{1}{2}$ " ball on holder



3 mm ball holder